

MR BREAST SOLUTIONS MADE BY NORAS

Biopsy Unit for BI 320:

Optional Components for Post&Pillar and Grid Biopsy System

Breast Biopsy Sets


Compatible with all NORAS BI 320 breast biopsy coils and patient rests

part number

Breast Biopsy Set Post&Pillar: autoclavable and reusable;
consisting of:

- Base plate 1 St.
- Base unit fixation plate 2 St.
- P&P positioning unit with telescope bar 1 St.
- P&P needle guide 1 St.
- P&P marker 1 St.
- Slat plate horizontal: medial and CC 2 St.
- Slat plate horizontal: lateral 1 St.
- Insertion plate 1 St.
- Manual BI320-PA 1 St.

118941



Breast Biopsy Set Grid disposable: autoclavable and reusable;
with disposable grids in sterile packaging*; **consisting of:**

- Basisplatte 1 St.
- Grundeinheit Fixierplatte 2 St.
- Grid Biopsie Einheit Höhenverstellbar: lateral (Einweg 5er Pack) 1 St.
- Grid Biopsie Einheit Höhenverstellbar: medial, cc (Einweg 5er Pack) 1 St.
- Grid Markerblock 1 St.
- Einsteckgitter horizontal - medial, cc 1 St.
- Einsteckgitter horizontal - lateral 1 St.
- Einlegeplatte 1 St.
- Gebrauchsanweisung BI320-PA 1 St.


118943



Breast Biopsy Set Grid reusable: autoclavable and reusable;
consisting of:









- Basisplatte 1 St.
- Grundeinheit Fixierplatte 1 St.
- Grid Biopsie Einheit Höhenverstellbar: lateral 1 St.
- Grid Biopsie Einheit: lateral, medial, cc 1 St.
- Grid Markerblock 1 St.
- Einlegeplatte 1 St.
- Gebrauchsanweisung BI320-PA 1 St.

118944









* distributed by NORAS

Optional Components for Post&Pillar Biopsy System

	part number	
Post&Pillar positioning unit with telescope bar: autoclavable and reusable	111259	
Post&Pillar Needle Guide: autoclavable and reusable; able to be angulated by 15° and 30°	111451	
Post&Pillar marker: autoclavable and reusable; (remove marker liquid before autoclaving and refill afterwards)	111300	
Post&Pillar needle guide sleeves*: disposable in sterile packaging <i>* distributed by NORAS</i>	18G 112654 16G 112597 14G 112655 12G 112653	 16G 18G 14G 12G
Base unit fixation plate: can be used laterally, medially and cranio-caudally and is used as adapter for slat plates and grids	111242	
Base plate	111292	
Slat plate horizontal: medial and cc	111212	
Slat plate horizontal: lateral	111301	

Optional Components for Grid Biopsy System

	part number	
Grid markerblock: autoclavable	111251	
Grid biopsy unit: lateral, medial, CC; autoclavable and reusable	111252	
Grid biopsy unit height adjustable: lateral; autoclavable and reusable; fits all BI320 systems (NORAS/Siemens)	117143	
Grid biopsy unit height adjustable*: lateral; disposable in sterile packaging; fits all BI320 systems (NORAS/Siemens); packaging consists of 5 grids <i>* distributed by NORAS</i>	112235	
Grid biopsy unit height adjustable*: medial, CC; disposable in sterile packaging; packaging consists of 5 grids <i>* distributed by NORAS</i>	112238	
Grid needle blocks for 18G, 16G, 14G and 12G*: disposable in sterile packaging <i>* distributed by NORAS</i>	<p>18G 112660</p> <p>16G 112143</p> <p>14G 112659</p> <p>12G 112731</p>	 <p>18G 16G</p> <p>14G 12G</p>

NORAS – Your partner for innovative solutions in MR Imaging

Since more than 30 years NORAS MRI products GmbH is a reliable partner for high resolution imaging in MR-based diagnostics and intervention. Kindly get in touch with us in order to learn more about the application areas of our solutions for your MR-based tasks.

NORAS MRI products GmbH
GERMANY
Leibnizstrasse 4
97204 Hoechberg
Tel.: +49 (0) 931 2 99 27 0
Fax: +49 (0) 931 2 99 27 20
info@noras.de
www.noras.de



NORAS MRI products GmbH
USA
Mr. John Wilkie
Tel.: +1 (414) 350 4433
sales-usa@mri-products.com
www.mri-products.com